

BLOOD PRESSURE MONITORING (0441)

1. Unstow BP Recorder Assembly
2. Replace Batt
3. RECORDER PWR - ON
4. Depress recessed RESET button on rear panel
5. Record MET
6. Set SAMPLE TIME on front panel to 3.5 min
7. Remove shirt
8. Rub electrode sites with tissue
9. Remove protective tape from (3) electrodes
10. Attach electrodes:

Red Lead - left lower chest
Green Lead - mid lower chest
White Lead - right lower chest

11. Remove protective tape from transducer adhesive 'donut'
12. Attach transducer to inner surface of non dominant arm 3" above elbow over brachial artery pulse
13. Wrap BP cuff around arm and transducer
14. Route leads, replace shirt
15. Plug in EKG, transducer cable and BP cuff hose
16. Attach recorder to belt, fasten around waist
17. Press MANUAL RECYCLE pb, collect two recording cycles
18. Verify system operating correctly

NOTE

To abort a cycle or initiate a new one, press MANUAL RECYCLE pb

19. When unit is removed:
RECORDER SAMPLE TIME - STBY

CAUTION
DO NOT TURN POWER OFF

MED

FLT

SUP

MED 3505

BOWEL SOUNDS RECORDING

1. Unstow Sony Walkman, cassettes and abdominal microphone belt
2. Remove shirt
3. Strap microphone belt around abdomen with Batt pack midline
4. Plug conector into Sony (red jack LINE IN/MIC)
5. Insert tape, record subject and MET on tape
6. Set recorder controls:
 - REC LEVEL - max clockwise
 - LINE IN/MIC - LINE IN
 - TAPE - NORMAL
 - DOLBY NR - OFF
 - COUNTER - 000
 - Simultaneously depress - PLAY, RECORD
7. Microphone batt pack - ON
8. Take 10mg. Reglan tablet, record MET
9. Continue recording until cassette complete (45 min)
10. If time allows reverse tape and continue for additional 45 min
11. Record MET and SUBJECT on second side of tape
12. Microphone batt - OFF
13. Disconnect and stow

TEMPERATURE, PULSE AND COLOR

1. Obtain paper thermometer from medical kit
2. Place under tongue at back of mouth
3. After 5 min remove, wait 1 min, record
4. Palpate radial artery at wrist, count 1 min, record
5. Obtain color chips
6. Move to source of good illumination (daylight if possible)
7. Place chips beside cheek and record nearest matched chip #

EKG RECORDING

1. Obtain Holter recorder
2. Record clock time in recorder window
 - flashing colon - PM
 - steady colon - AM
3. Grasp silver electrode handles (pb in right hand)
4. Depress pb for 1 min
5. Stow recorder

MED
FELT / MED DESO

MCP ACTIVITY DATA SHEET

SUBJECT					
MET REGLAN TAKEN					
TEMP					
PULSE					
COLOR					

SUBJECT					
MET REGLAN TAKEN					
TEMP					
PULSE					
COLOR					

SUBJECT					
MET REGLAN TAKEN					
TEMP					
PULSE					
COLOR					

SUBJECT					
MET REGLAN TAKEN					
TEMP					
PULSE					
COLOR					

MED

100 0505